



Naziv uređaja Pumpa za transkardijalnu perfuziju	Apparatus Transcardial perfusion pump
Proizvođač i model/The manufacturer and model	
Kratak opis metode Pumpa imitira rad srca i postepeno uvodi perfuzioni rastvor u srce životinje koja je pod anestezijom. Brzinu uvođenja rastvora je moguće podešavati.	Short description of the method The pump is imitating heartbeats and it gradually infuses a perfusion solution into the heart of the anaesthetized animal. The infusion speed is adjustable.
Tehničke karakteristike Elastična creva malog promera povezuju pumpu sa perfuzionom iglom i sudom u kojem je perfuzioni rastvor. Pumpa sadrži prekidač pomoću koga je moguće podesiti brzinu protoka iz suda iz koga se uzima perfuzioni rastvor.	Technical characteristics Elastic tubes connect the pump with the perfusion needle and the dish with perfusion solution. It is possible to set the flow speed with the switch on the pump.
Primena i tip uzorka Imobilisana eksperimentalna životinja koja je pod anestezijom	Application and sample type Immobilized anaesthetized experimental animal
Osoba za kontakt / Contact person Doc. dr Marin Jukić, marin.jukic@pharmacy.bg.ac.rs	
Link ka uređaju na sajtu proizvođača / Link of the product on the manufacturer's website	