



<b>Naziv uređaja</b> Vakuum koncentrator RVC 2-18 CDplus	<b>Apparatus</b> Vacuum concentrator RVC 2-18 CDplus
<b>Proizvođač i model/The manufacturer and model</b> Martin Christ (Germany), RVC 2-18 CDplus	
<b>Kratak opis metode</b> Uparavanje uz rotaciju	<b>Short description of the method</b> Evaporation with rotation
<b>Tehničke karakteristike</b> Brzina rotora do 1500 rpm, Vacuubrand membranska pumpa MZ 2C (postignuti vakuum 8 mbara), 2 raspoloživa rotora (veliki: za zapremine od 10 i 15 mL, kapacitet max 12 tuba; mali: za zapremine od 1,5 mL, kapacitet max 24 Eppendorf tube)	<b>Technical characteristics</b> Speed rotor 1500 rpm, Vacuum pump MZ 2C (end vacuum 8 mbar), 2 available rotors (large: for volumes of 10 and 15 mL, capacity max 12 tubes; small: for volumes of 1,5 mL, cašcapacity max 24 Eppendorf tubes)
<b>Primena i tip uzorka</b> Delimično ili uparavanje do suva organskih i vodenih rastvora	<b>Application and sample type</b> Partial or evaporation to dry organic and aqueous solutions
<b>Osoba za kontakt / Contact person</b> Tamara Antičić	
<b>Link ka uređaju na sajtu proizvođača / Link of the product on the manufacturer's website</b> <a href="https://www.martinchrist.de/en/products/rvc/rvc-2-18-cdplus/">https://www.martinchrist.de/en/products/rvc/rvc-2-18-cdplus/</a>	